



Faux Leather/Vinyl Specifications

RANGE	Seawater Resistant	Chlorine water resistant	Resistant to sweat	Resistant to urine & blood	OKO-Tex Standard 100	EN 71-3 safe toys	Phthalate free	DINEN ISO 10993 (Part 5+10) medical product law	Antibacterial properties	DIN EN ISO 105-B02 light fastness	Abrasion (martindale)	Thickness in mm	Weight in g/m ²	Width in cm	Length in mm
Basis										6	50,000	1.3 ± 0.1	820 ± 50	140	25/30
Baroness			✓				✓			6	100,000	1.2 ± 0.1	610 ± 30	140	30
Countess	✓	✓	✓	✓	✓	✓	✓			6	75,000	1.1 ± 0.1	650 ± 40	140	25
Lady	✓				✓					6	100,000	1.1 ± 0.1	610 ± 50	140	25
Marquis					✓					6	100,000	1.0 ± 0.1	540 ± 50 w	140	25
Prince										5	305,000	1.45 ± 0.2	535 ± 50	140	30
Princess				✓			✓			6/7	50,000	1.2	610 ± 30	140	30
Queen			✓	✓	✓	✓	✓	✓		6	50,000	1.1 ± 0.1	650 ± 40	140	25
Viscount	✓	✓	✓	✓	✓	✓	✓		✓	6	100,000	1.2 ± 0.1	740 ± 50	140	25

Please ask for datasheets for tested disinfectant agents.

Processing: Trouble free handling in usual upholstery industry processes (sewing, tacking, gluing, nailing, HF-welding). No solvent adhesives to be used.

Cleaning: Soiled areas can be cleaned with a damp cloth. Particularly recommended is the use of a commercially available microfibre cloth. Stubborn stains require treatment with a commercial mild action detergent.

What you should know about ordering faux leather

Samples are a good guide to colour and finish, but you should bear in mind an exact match cannot be guaranteed.

Email | sales@ukhide.co.uk

Phone | +44 (0)1932 353 338

Visit us online at

www.ukhide.co.uk